

WHAT WE NEED TO MADE A PERFECT NEW SMILE FOR YOUR PATIENT

Before unscrewing the conversion:

- 1 Take photos of:
 - a Midline
 - b Entire face, smiling
 - c Retracted view pt biting completely closed
 - d Profile photos while smiling AND retracted
 - e Custom shade tab against the patients teeth or converted denture
 - f Short video of the lip dynamics while smiling and talking
- 2 Easy Bar Impression using titanium cylinders. YouTube video:
<https://www.youtube.com/watch?v=1fbf2OrDZFA>
- 3 Open Tray Impression of Easy Bar capturing the intaglio surface
- 4 Duplication of the existing converted denture YouTube video showing how:
<https://www.youtube.com/watch?v=H4CEoO6ladj&list=PLFygAIxWTixJvq9k2lFCGiL8R-pm2i->
- 5 Bite Registration
- 6 One set of lab analogs AND titanium cylinders (depending on the package that you choose)

For the economy package we will need one set of cylinders and if the patient wants to WEAR the PMMA, we will need 2 sets of non engaging titanium cylinders